




I am Armadillo.

I'm a new sort of catheter for a transradial approach. My radial-access system is elegantly tailored to the demands of accessing neurovasculature through the radial pathway. Plus, I'm optimized with two modes for seamless adaptation.



On the path to radial-access solutions, I'm leading the way.

The radial access approach is revolutionizing what's possible in neurovascular procedures—and I'm a front-line partner for caregivers propelling the profession forward. With my dual-mode capabilities and slick construction, I'm a solid option for the patient and a flexible one for the physician.



Tip Shape Straight	Outer Diameter (OD) 0.095"	Guidewire Compatibility 0.035" - 0.038"
------------------------------	--------------------------------------	---

0.072" ID

My 0.072-inch internal diameter enables procedural flexibility.

95cm, 105cm

My length and SelectFlex technology mean I'm a game-changer when it comes to trackability and support.


7 French

My 0.095-inch OD allows for compatibility with a 7 French system.

Accessories

- Mode Control Syringe
- Peel-Away Introducer
- Luer-Activated Valve
- Dilator

Here's how I deliver.



HOW I'M DIFFERENT

I'm nimble but solid.

When it comes to radial access, a slim profile is key—but so is superior support. With my kink-resistant construction, powerful coil design, and SelectFlex variable stiffness technology, I can navigate even the most delicate tortuosity, providing a tremendously stable platform.

Attributes:

- Dual-mode technology with tracking and support modes
- 22- or 27-centimeter variable support transition length options


HOW I'M DIFFERENT

I'm a smooth operator.

With an all-in-one system and accompanying dilator accessory, I navigate acute bends with strength that outstrips my size and agility in the slimmest spots. How do I do it? With transition zones generated specifically for navigating through the anatomy from a radial approach.

Attributes:

- Elastically loaded tapered tip for a zero gap transition
- Dilator lengths vary with catheter length and have a luer-fitting connection



HOW I'M DIFFERENT

I'm actually quite advanced.

Fluidity is fundamental when it comes to navigating some of the most delicate blood vessels in the body. So I'm armed with a hydrophilic coating for optimal movement without resistance.

Attributes:

- 30cm of hydrophilic coating
- Lubricious and highly durable

Product Code	Description	Length	Variable Stiffness Length
SF7072-095-22-30	Selectflex, Armadillo 7F	95cm	22cm
SF7072-095-27-30	Selectflex, Armadillo 7F	95cm	27cm
SF7072-105-22-30	Selectflex, Armadillo 7F	105cm	22cm
SF7072-105-27-30	Selectflex, Armadillo 7F	105cm	27cm

We're Q'Apel, a medical technology company that creates solutions. More specifically, we design novel access device technology for vascular interventions and unmet clinical needs. Because in the precious seconds that surround a stroke emergency, clinicians need technology that delivers. That's where we come in.

Get in touch | info@qapelmedical.com | [linkedin.com/company/q-apel-medical](https://www.linkedin.com/company/q-apel-medical)

©2021. Q'Apel Medical, Inc. All rights reserved. MKT 00621.D

Q'Apel Medical
46708 Lakeview Blvd
Fremont, CA 94538

CALL TO INQUIRE: 510.738.6255

[qapelmedical.com](https://www.qapelmedical.com)